

# RePo-1, cat. no. 5010



## Portable Refracto-Polarimeter

### Fusion of refractive index and angle of optical rotation

Brix alone can only indicate the overall amount of dissolved solids in a solution. With the RePo-1 it is now possible to know the actual quantities of each dissolved solid, such as fructose and glucose content, in a solution.

## RePo Four Points

**Point 1 – 2 in 1** - Brix alone can only indicate the overall amount of dissolved solids in a solution. With RePo-1, it is possible to determine if a product has been adulterated by measuring the optical rotation. For example, fructose rotates light to the left while glucose rotates light to the right. If a sample has been adulterated with glucose, it will make the optical rotation readings appear more positive than a pure solution.

**Point 2 - Compact size** - Great portability, 10.1 × 16.0 × 3.8cm and weighs only 325g.

**Point 3 - Automatically calculated scales** - Press "User" to select the desired scale: Purity, International Sugar Scale, Specific Rotation, and Concentration.

A.R. –Angle of Rotation



Brix - Brix



User Specific Scales



**Point 4 - Measurement evaluation indicator** - Measurements are evaluated against pre-set tolerance limits, and the result is indicated by green or red light. (Upper and lower tolerance limits must be set in advance to using this function)



Green light - Within tolerance



Red light - Outside tolerance

## RePo-1 Application Examples

- ❖ Sugar industries – sugar manufacturers and refineries
- ❖ Analyzing inverted sugar syrup
- ❖ Analyzing fragrances
- ❖ Honey
  - ❖ Investigating how fast different honey crystallizes by determining the fructose and glucose content
  - ❖ Researching honey collected from the same district or honey collected from the same flower

## Measurement Method

1. Fill the sample stage up to the fill line (approx. 3mL)



2. Press START.



3. Measurement is displayed



## RePo-1 Specifications

Scale	Angle of rotation, Brix, International Sugar Scale, Temperature
Measurement Readings	Angle of rotation, Temperature, User (International Sugar Scale/ Purity/ Specific Rotation/ Concentration), Brix, Purity, International Sugar Scale with ATC
Measurement Range	Angle of Rotation: -5.00 to +5.00°(*1) International Sugar Scale: -130 to +130°Z Brix: 0 to 85% Temperature: 15.0 to 40.0°C
Display Range	Angle of Rotation: -5.99 to +5.99° Brix: -2.0 to 86.6% International Sugar Scale: -130 to +130°Z Purity: -120 to 120% Specific Rotation: -999 to +999 Concentration: -2.0 to 999% Temperature: 14.0 to 41.0°C
Resolution	Angle of Rotation: 0.01° Brix: 0.1% International Sugar Scale: 0.1°Z Temperature: 0.1°C
Measurement Accuracy	Angle of Rotation: ±0.1°(at 20°C) International Sugar Scale: ±3°Z Brix: ±0.2% Temperature: ±1°C
Repeatability	Angle of Rotation: ±0.05°(at 20°C) Brix : ±0.1% International Sugar Scale: ±1.5°Z (at 20°C)
Wavelength	589nm (equivalent to the sodium-D line spectrum)
Temperature compensation range	Brix: 15 to 40°C Purity: 18 to 40°C International Sugar Scale: 18 to 40°C

### RePo Series - In addition to angle of rotation and Brix:

**RePo-2:** can determine the amount of fructose (%) in HFCS - ideal for various industries such as food and beverage, or sugar refining.

**RePo-3:** can determine the amount of invert sugar (%) converted from sucrose - ideal for various industries such as sugar refining.

**RePo-4:** can determine the amount of fructose (%) and honey moisture (%) - ideal for Q.C. of honey.

**RePo-5:** can determine refractive index (nD) - ideal for quality control in various industries such as fragrances, cosmetics, or pharmaceutical.

t: 01580 891300

e: qa@aceindustrial.co.uk